



IV Catheter

Intravenous Ethylene Tetrafluoroethylene (ETFE)



- Designed to perform in delicate vasculature, delivering maximum efficiency and patient comfort
- Smooth catheter-to-needle transition results in easy insertion
- ETFE catheter material is biocompatible, minimizing risk of irritation and infection while catheter is indwelling
- Hard plastic outer packaging protects catheter before use
- Radiopaque catheter features flashback “window” for quick visualization of blood
- Hub is color-coded for rapid identification of catheter size

ORDERING INFORMATION

| Item # | Description | Packaging |
|---------|-------------|--------------------|
| CI+1450 | 14GX2 | |
| CI+1632 | 16GX1.1/4 | |
| CI+1650 | 16GX2 | 50/Box |
| CI+1832 | 18GX1.1/4 | 4 Boxes/Case (200) |
| CI+1850 | 18GX2 | |
| CI+2025 | 20GX1 | |

| Item # | Description | Packaging |
|---------|-------------|--------------------|
| CI+2032 | 20GX1.1/4 | |
| CI+2050 | 20GX2 | |
| CI+2225 | 22GX1 | 50/Box |
| CI+2232 | 22GX1.1/4 | 4 Boxes/Case (200) |
| CI+2420 | 24GX3/4 | |
| CI+2425 | 24GX1 | |

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



NIPRO MEDICAL CORPORATION
200 Crossing Blvd., Bridgewater, NJ 08807
T: +1 (908) 393-7030, F: +1 (908) 393-7031
www.nipro.com

© 2014 Nipro Medical Corp. All Rights Reserved.
Printed in U.S.A. MS-DSUS-IVCATH Rev 0