




NiproSet[®] Blood Tubing



- Built-in safety and efficiency features, with all sets pre-labeled based on equipment compatibility


- Demonstrated accurate and consistent flow rates as well as excellent performance under extreme conditions, when compared with a competitive blood tubing set



NIPRO CONNECT[®] PROGRAM

Our commitment to strengthen every partnership.

All Nipro renal products are backed by the Nipro Connect Educational Program – offering clinicians product evaluations, in-services, continuing education units (CEU's) and instructional materials.



Medication ports

Large finger guard helps protect the user from accidental needle sticks; venous line features multiple medication ports for optimum versatility and choice in medication delivery.



Color-coded clamps

Patient safety is enhanced with clamps that are color-coded for confident identification of the appropriate line.




Transducer protectors

Factory-applied torque helps prevent loosening or disconnecting of transducer; engineered to maintain monitoring while decreasing exposure to the dialysis machine.



Drip chambers

Drip chambers feature a versatile medication port and, for an extra measure of safety, are anti-foaming to minimize splashing.



In-vitro testing shows that when compared with a leading brand of blood tubing, NiproSet delivered an accurate and consistent flow rate over time.*

*Data on file

ORDERING INFORMATION

Item #	Description	Packaging
BL+A204/V801	Baxter1000/Baxter Arena/ 8mm/Post-Pump	24/Case (24)
BL+A209Y/V803	Fresenius 2008 series/ 8mm/Pre-Pump	
BL+A210/V803	Fresenius 2008 series/8mm/ Post-Pump	
BL+A211/V804	Baxter 550-1550/8mm/ Post-Pump	

Item #	Description	Packaging
BL+A214Y/V805	Fresenius 2008 series/8mm/ No Arterial Chamber	24/Case (24)
BL+A217Y/V806	B. Braun Dialog/8mm/ Post-Pump	
BL+A223Y/V809	Fresenius 2008 series/8mm/ Pre-Pump/Pre-Attached Priming Set	

CAUTION: Federal law restricts these devices to sale by or on the order of a physician.



NIPRO MEDICAL CORPORATION
200 Crossing Blvd., Bridgewater, NJ 08807
T: +1 (908) 393-7030, F: +1 (908) 393-7031
www.nipro.com

© 2014 Nipro Medical Corp. All Rights Reserved.
NiproSet is a registered trademark of Nipro Corporation.
Printed in U.S.A. RN-DSUS-BLTUB Rev 0