



# Infusion and Blood Collection Systems



All Nipro Infusion and Blood Collection Systems are designed to deliver three important benefits:


- Maximize user safety
- Minimize patient discomfort
- Enhance ease of use

Blood Infusion and Collection products:


- SafeTouch® Winged Infusion Sets
- Winged Infusion Sets
- Blood Collection Needles
- Luer Adapter

# Blood Infusion

## SafeTouch Winged Infusion Sets



## Winged Infusion Sets



- The wings are soft and color-coded with a non-slip grip.
- Patient comfort is assured with a needle wall that is significantly thinner than a regular wall and virtually painless.
- Infusion set tubing is flexible.
- SafeTouch has built-in safety features that help minimize the incidence of accidental needlesticks.

## ORDERING INFORMATION

### Nipro SafeTouch Winged Infusion Set

Item #	Description	Packaging
<b>SPR+19G19</b>	19GX3/4 12" Tubing	50/Box, 10 Boxes/Case (500)
<b>SPR+21G19</b>	21GX3/4 12" Tubing	
<b>SPR+23G19</b>	23GX3/4 12" Tubing	
<b>SPR+25G19</b>	25GX3/4 12" Tubing	

### Nipro Winged Infusion Set

Item #	Description	Packaging
<b>PR+18G19</b>	18GX3/4 12" Tubing	50/Box, 20 Boxes/Case (1,000)
<b>PR+19G19</b>	19GX3/4 12" Tubing	
<b>PR+21G19</b>	21GX3/4 12" Tubing	
<b>PR+22G19</b>	22GX3/4 12" Tubing	
<b>PR+23G19</b>	23GX3/4 12" Tubing	
<b>PR+25G19</b>	25GX3/4 12" Tubing	
<b>PR+27G19</b>	27GX3/4 12" Tubing	

**CAUTION:** Federal law restricts these devices to sale by or on the order of a physician.

# Blood Collection

## Blood Collection Needle



- Multiple-Draw needle design provides the convenience of several sample collections with a single venipuncture.
- Needle is ultra-sharp and with smooth edges for ease of penetration, resulting in atraumatic insertion.
- Color-coded hubs assure accurate and convenient gauge size identification.

## Luer Adapter



- Allows for multiple draws with a single venipuncture.
- Designed to accommodate a wide variety of devices.

## ORDERING INFORMATION

### Nipro Blood Collection Needle

Item #	Description	Packaging
NM+20G25	20GX1	
NM+20G38	20GX1.1/2	
NM+21G25	21GX1	100/Box, 10 Boxes/Case
NM+21G38	21GX1.1/2	(1,000)
NM+22G25	22GX1	
NM+22G38	22GX1.1/2	

### Nipro Luer Adapter

Item #	Description	Packaging
LM+20G	Multi Sample Luer Adapter	100/Box, 10 Boxes/Case (1,000)

**CAUTION:** Federal law restricts these devices to sale by or on the order of a physician.



**NIPRO MEDICAL CORPORATION**  
200 Crossing Blvd., Bridgewater, NJ 08807  
T: +1 (908) 393-7030, F: +1 (908) 393-7031  
[www.nipro.com](http://www.nipro.com)

© 2014 Nipro Medical Corp. All Rights Reserved.  
SafeTouch is a registered trademark of Nipro Medical Corporation.  
Printed in U.S.A. MS-BRUS-INF Rev 0